

# Western Blotting

92/5/20

written by Tomoaki Shirao ver 1.03

## Blotting buffer (Laemmli - SDS)

2 x Tris-glycine (ph8.5)      500 ml  
MeOH                              200 ml  
add water to 1000 ml

## Procedure

1. remove the gel from the plate
2. immerse it into the blotting buffer
3. sandwich the nitrocellulose sheet and the gel between filter sheets (Those sheets should be preincubated.)
4. Put the above sheets in the blotting apparatus: red side is nitrocellulose side (red is plus)
5. 200 mA, overnight for drebrin transfer

## Staining with amidoblack

staining solution:      45% MeOH                      45 ml  
                                 10% AcOH                      10 ml  
                                 0.1% Amido Black              0.1g  
                                 add water to                      100 ml

## Immunostaining

1. incubation in 10% milk                      60 min
2. brief PBs Wash
3. first antibody in PBSA (3% albumin)      60 min
4. PBS wash                      5 min
5. incubation in 10% milk                      5 min
6. brief PBS wash
7. second antibody in in PBSA (3% albumin)      60 min  
(20 µl in 5 ml)
8. PBS wash                      15 min x 3
9. reaction

## HRP reaction

solution A:                      0.3% chloronaphthol (12 mg in 4 ml MeOH)  
solution B:                      0.3% hydrogen peroxide (x 1/100)  
solution C:                      imidazole 27 mg/ml

## Procedure

1. 20 ml of PBS + 4 ml of A + 0.8 ml of B + 0.5 ml of C
2. wash in water